



PhoenixNMR Does DNP!

Introducing the PhoenixNMR 2.5 mm 18 μ l sample volume 395 GHz / 600 MHz $^1\text{H}/^{19}\text{F}/(^{31}\text{P}-^{13}\text{C})/(^{23}\text{Na}-^{15}\text{N})$ quadruple resonance MAS-DNP NMR probe. Microwave power is transported through probe using a series of corrugated miter and straight sections connected to a focusing lens designed by Thorsten Maly³ of Bridge12 Technologies. At a modest cost well below any other complete DNP probe system on the market, with greatly expanded operational capability including a low gamma box to reach all NMR accessible nuclei, this probe is a perfect addition to your existing DNP system (400-600 MHz for NMR) or simply add a Bridge12 gyrotron to an existing NMR system to have a full DNP system.





DNP Enhancement versus sample temperature: The four round markers reflect the data collected on the PhoenixNMR 2.5 mm DNP probe utilizing the 395 GHz/600MHz DNP NMR system at the NHMFL¹, VT flow rate of 30 lpm, the two orange dots represent the performance of the Bruker 3.2 mm DNP probe on the same system with the same sample, the single green dots shows the expected performance based on recent improvements to the VT system of the PhoenixNMR probe at 40 lpm flow rate.

^{13}C - ^{13}C 2D correlation of ^{13}C labeled intact plant stem (Sample and Data by Tuo Wang²) on the 600 MHz DNP system at the National High Magnetic Field Laboratory in Tallahassee, FL, in collaboration with Frederic Mentink-Vigier¹. This spectrum was acquired with a spin rate of 22 kHz, a sample temperature of $\sim 140\text{K}$ (run at a conservative 20 lpm VT flow) exhibiting a DNP enhancement of 8:1.



Sample temperature measurement at a set temperature of 90 K utilizing the chemical shift of PbNO_3 . The central 5 spectra show the sample temperature as a function of spin rate at a VT flow rate of 30 lpm. The top and bottom most spectra show the effect of increasing the VT Flow from 30 lpm to 40 lpm.






The PhoenixNMR DNP probe system: The WB Quad Resonance HFX probe (a), two control boxes (b,c) - MAS spinning, VT cryogenic control and sample insert/eject operation), cryogenic counter-flow heat exchanger (d), LN₂ Dewar (not Shown). The 50-liter Dewar provides ~36 hours of 90 K operation and can be refilled at any time. The included counter-flow heat exchanger is manufactured by Bridge12 Technologies³.

10mM AMUPOL, DNP Juice, 0.25M ¹³C-¹⁵N Proline
 18μL sample, T=150K, μWave ~ 10 W
 Enhancement: 42:1
 MAS = 25 kHz

Blue = μWave on
 Red = μWave off



Tuning Range

Quadruple Resonance Mode: **H:** ¹H and ¹⁹F **X:** ³¹P to ¹³C **Y:** ²³Na to ¹⁵N
 Triple Resonance Mode: **H:** ¹H and ¹⁹F **X:** ³¹P to ¹⁵N
 Tuning on X can be lowered to ~15 MHz with an optional Low Gamma box.
 Simultaneous H and F operation available as an option.

Spinning Module

Spin Rate: 5 kHz to 25 kHz @ 90 K to 373 K (set temperature)
 Stability: +/- 5 Hz or +/- 0.2%, whichever is larger, over 24 hrs
 Sample Volume: 18 μL Active sample volume used for all NMR tests
 VT Range: 90 K to 373 K
 90 K @ 5 kHz Spin Rate, Sample T = 115 K
 90 K @ 25 kHz Spin Rate, Sample T = 130 K

Resolution

Adamantane (¹³C): 0.05 ppm FWHM, 0.2 ppm @ 10%, 0.35 ppm @ 2%

Signal to Noise (Glycine, Matched filter, 32 acquisitions, 293 K)

HFCN Mode 150:1
 HFC Mode 170:1

RF Performance @ 90 K and 300 K:

Nucleus	pw90μs
¹ H/ ¹⁹ F	2.5 μs (50 ms, 2% duty factor)
¹ H & ¹⁹ F	2.8 μs (50 ms, 2% duty factor)
¹³ C	2.5 μs (10 ms, 2% duty factor)
³¹ P	2.5 μs (10 ms, 2% duty factor)
¹⁵ N	5.0 μs (10 ms, 2% duty factor)

¹ Frederic Mentink-Vigier, *National High Magnetic Field Laboratory, Tallahassee, Florida 32310*, ² Tuo Wang, *Department of Chemistry, Louisiana State University, Baton Rouge, Louisiana 70803*, ³ Thorsten Maly, *Bridge12 Technologies, Inc. Framingham Massachusetts 01702*

The National High Magnetic Field laboratory (NHMFL) is funded by the National Science Foundation Division of Materials Research (DMR-1157490 and 1644779) and the State of Florida. The DNP system at the NHMFL is funded in part by NIH S10 OD018519, P41 GM122698, and NSF CHE-S1229170.