



PhoenixNMR Does DNP!

Dynamic Nuclear Polarization (DNP) signal-enhancement capability is now available in combination with innovative Phoenix Solid-State NMR (SSNMR) probe technology to expand the scope and sensitivity of biomolecular, pharmaceutical, and materials SSNMR research at 400-600 MHz. The new Phoenix MAS-DNP probe system utilizes the same robust, professional design strategy that scientists have come to expect in all Phoenix SSNMR probes and accessories. Bridge12 Technologies components pair with Phoenix probes and controllers to provide the full DNP system.



The PhoenixNMR MAS-DNP probe system: a) WB Quad Resonance HFX probe, b) MAS controller, c) VT cryogenic and sample insert/eject controller, d) Cryogenic counter-flow heat exchanger. Not shown: 50-liter LN₂ Dewar. The 50-liter Dewar provides ~36 hrs of operation at T = 90 K and can be refilled at any time. The included counter-flow heat exchanger is manufactured by Bridge12 Technologies³.



¹³C CP/MAS DNP enhancement of 42:1 was measured^{1,2} using a PhoenixNMR 395 GHz/600 MHz DNP HFX probe fitted with a 2.5 mm spinning module. Sample Temperature = 150 K at MASS = 25 kHz. The sample was prepared using 10 mM AMUPol biradical and "DNP Juice" = d₈-Glycerol/H₂O/D₂O in a 6/3/1 ratio.

Overview of the Phoenix DNP Innovation

- Add a new Phoenix DNP probe to expand capability on an existing DNP SSNMR instrument; or
- Add a Turn-key DNP System, including gyrotron, to an existing 9.4 – 14.1 T, Wide-Bore SSNMR instrument, pairing the robust Phoenix NMR probe technology with Bridge12 Technologies components to complete the DNP package;
- Combine DNP with HFX and HFX probe configurations;
- Optional Probe Accessories allow low-gamma nuclei and/or simultaneous H/F operation;
- Operating Temperature Range is 90 K to 373 K (actual sample temperature is dependent on spinning speed);
- 2.5 mm MASS Module is available for spinning speeds up to 25 kHz;
- Insert and Eject Samples In-Situ at cold and warm temperatures;
- Modest up-front Costs beat other commercially-available DNP probe systems;
- Unique VT Design minimizes the use of cryogenics to dramatically reduce operational costs.

Full DNP Probe System from Phoenix NMR includes:

- Phoenix NMR DNP wide bore (WB) HFX or HFX probe (low-gamma accessory, simultaneous H and F operation, and lock channel optional) including sample eject;
- Phoenix NMR MAS speed controller;
- Phoenix NMR VT cryogenic and sample insert/eject controller;
- 50-liter LN2 dewar, provides ~36 hours of operation per fill at 90 K, and can be refilled during operation;
- Bridge12 Technologies Waveguide utilizing corrugated waveguide sections and a Focusing Lens directed perpendicular to sample coil for maximum DNP efficiency;
- Bridge12 Technologies Cryogenic counter-flow heat exchanger.

Full DNP System from Phoenix NMR includes the above, plus:

- Bridge12 Technologies fastwave gyrotron operating at 263 GHz (for 400 MHz instrument) or 395 GHz (for 600 MHz instrument) and controller.

PhoenixNMR DNP Probe Performance*

Tuning Range

Quadruple Resonance Mode: **H:** ^1H and ^{19}F **X:** ^{31}P to ^{13}C **Y:** ^{23}Na to ^{15}N

Triple Resonance Mode: **H:** ^1H and ^{19}F **X:** ^{31}P to ^{15}N

Tuning on X can be lowered to ~15 MHz with optional Low Gamma box.

H channel can be tuned to simultaneously deliver $^1\text{H}/^{19}\text{F}$ power and/or observe $^1\text{H}/^{19}\text{F}$ with an optional probe accessory.

Spinning Module

Spin Rate: 5 kHz to 25kHz @ 90 K to 373 K (set temperature)

Stability: +/- 5 Hz or +/- 0.2%, whichever is larger, over 24 hrs

Sample Volume: 18 μL Active sample volume used for all NMR tests

VT Range: 90 K to 373 K

90 K @ 5 kHz Spin Rate, Sample T = 115 K

90 K @ 25 kHz Spin Rate, Sample T = 130 K

Resolution

Adamantane (^{13}C): 0.05 ppm FWHM, 0.2 ppm @ 10%, 0.35 ppm @ 2%

Signal to Noise (Glycine, Matched filter, 32 acquisitions, 293 K)

HFCN Mode 150:1

HFC Mode 170:1

RF Performance @ 90K and 300K:

<u>Nucleus</u>	<u>pw90μs</u>
$^1\text{H}/^{19}\text{F}$	2.5 μs (50 ms, 2% duty factor)
^1H & ^{19}F	2.8 μs (50 ms, 2% duty factor)
^{13}C	2.5 μs (10 ms, 2% duty factor)
^{31}P	2.5 μs (10 ms, 2% duty factor)
^{15}N	5.0 μs (10 ms, 2% duty factor)

*Probe Performance Specifications are quoted herein for a Phoenix NMR 395 GHz/600 MHz DNP HFX probe fitted with a 2.5 mm spinning module, utilizing a Bridge12 Technologies gyrotron, waveguide, and counter-flow heat exchanger. Other hardware configurations may result in different probe performance specifications. Request a quote for detailed specifications.

¹ Frederic Mentink-Vigier, *National High Magnetic Field Laboratory, Tallahassee, Florida 32310*, ² Tuo Wang, *Department of Chemistry, Louisiana State University, Baton Rouge, Louisiana 70803*, ³ Thorsten Maly, *Bridge12 Technologies, Inc. Framingham Massachusetts 01702*

The National High Magnetic Field laboratory (NHMFL) is funded by the National Science Foundation Division of Materials Research (DMR-1157490 and 1644779) and the State of Florida. The DNP system at the NHMFL is funded in part by NIH S10 OD018519, P41 GM122698, and NSF CHE-S1229170.